

MSM® DEPERI

Endovascular PTA Balloon Catheter OTW 0.035" GW - Paclitaxel coatet



New tip and balloon

Advantages of the balloon

New balloon folding with lower profile

- excellent refolding, controlled compliance, rapid deflation
- Balloon diameter from 4-12mm
- Balloon length form 20-200mm



Large portfolio up to 200mm

Advantage of the catheter system

Dual Lumen shaft technology

- High kinking resistance
- Good track and pushability
- Easy maneuverability in tortous vessels