


MSM® DLT Deep

RX 2.8F Distal Limb Treatment 0.014"GW



2.8F shaft design




easy navigation

Benefits of the Balloon Catheter

- 2.8F shaft / Monorail System
- Balloon length from 10 to 40mm
- Balloon diameter from 1.5 to 4.0mm
- Excellent pushability
- High kink resistance
- Hydrophilic coated
- Very small lesion entry and crossing profile

Advanced Coating Technique

- durable hydrophilic coating on entire balloon
- no crack during deployment



Technical Details Catheter

Usable length		138 cm
Length of the guide wire lumen		27 cm
Balloon material	Balloon	Polyamide-Compound
	Distal shaft	Polyamide, multi-layer tubing
	Proximal shaft	Stainless Steel PTFE coating
Shaft size	Proximal	1.9 F
	Distal	2.8 F
Balloon folding	all sizes	3-fold
X-ray marker		Platinum / Iridium
Tip profile		0.017" (0.43 mm)
Max. recommended guide wire		0.014" (0.36 mm)

Size and Order Matrix

Balloon Length					
Balloon Ø	10 mm	15 mm	20 mm	30 mm	40 mm
1.5 mm	DLT 15010	DLT 15015	DLT 15020	DLT 15030	-
2.0 mm	DLT 20010	DLT 20015	DLT 20020	DLT 20030	DLT 20040
2.25 mm	DLT 22510	DLT 22515	DLT 22520	DLT 22530	DLT 22540
2.5 mm	DLT 25010	DLT 25015	DLT 25020	DLT 25030	DLT 25040
2.75 mm	DLT 27510	DLT 27515	DLT 27520	DLT 27530	DLT 27540
3.00 mm	DLT 30010	DLT 30015	DLT 30020	DLT 30030	DLT 30040
3.5 mm	DLT 35010	DLT 35015	DLT 35020	DLT 35030	DLT 35040
4.00 mm	DLT 40010	DLT 40015	DLT 40020	DLT 40030	DLT 40040